FILLMED

LABORATOIRES

NANOSOFT®
MICRONEEDLES

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

**DEVICE PREPARATION
FOR INJECTION**

- Draw injection solution into the syringe without purging.
- Open the soft blister and **firmly connect** the NANOSOFT® to the syringe.
- Remove the plastic cover.
- Purge air and excess fluid leaving the required volume in the syringe.
- Make sure injection site is dry prior to injection.

- **Hold syringe flanges** between thumb and fingers.
- The **blue line** must face you at all times.
- Syringe flanges are not the orientation guide as they may not always be parallel to the blue line; use the blue line for orientation.

- **Stretch the skin** downwards throughout the entire injection.
- Insert the device at a **45° angle**.
- Press the device into the skin until a **shallow indentation** is formed.

- **Inject slowly**, expect strong resistance.
- **A white wheal (bleb)** will form at the injection site.

NANOSOFT®

Microneedles for intradermal injection

For any further assistance/support/questions, please contact the manufacturer.

**MANUFACTURER**

NanoPass Technologies Ltd.

3 Golda Meir St., Nes Ziona 7403648, Israel

Tel +972-8-9462905

E-mail: info@nanopass.com

EC REP**EU REPRESENTATIVE**

MedNet EC-REP GmbH

Borkstrasse 10-48163, Muenster, Germany

Phone 49-251-322660

Fax 49-251-3226622

EN

NANOSOFT® MICRONEEDLES

NANOSOFT® INSTRUCTIONS FOR USE

DESCRIPTION

The NANOSOFT® is a microneedle device that is mounted on a standard syringe in the same manner as a conventional needle. The NANOSOFT® comprises an integral array of 3 microneedles; each microneedle 0.6mm in length. The NANOSOFT® is used to inject liquid substances, allowing for controlled intradermal delivery in any procedure which requires administration of substances to the dermal compartment. The NANOSOFT® is latex free.

INDICATIONS

The NANOSOFT® is intended for intradermal injections of any substance or drug approved for delivery by this delivery route.

CONTRAINDICATIONS AND LIMITATIONS

- The NANOSOFT® should not be used on skin abrasions, open wounds, cuts and scars.
- The NANOSOFT® should not be used on rashes, skin infections or any other area of damaged or diseased skin.

WARNINGS AND PRECAUTIONS

- The NANOSOFT® is intended for single use, by healthcare professionals only.
- The NANOSOFT® can be used in all ages and populations.
- Re-use may lead to infection or other illness/injury.
- Used or un-used devices should not be re-capped or removed from syringe unless there is no alternative or such action is required by a specific medical procedure.
- The NANOSOFT® is not intended for aspiration of liquids.
- If the device appears to be damaged or broken, discard and use a new device.
- Any serious incidence should be reported to the manufacturer and the competent authority.

ADVERSE EVENTS

Local adverse reactions, including edema, erythema and discoloration of skin at the injection site have been observed, most likely related to the injected substance, and not to the NANOSOFT®. Adverse reactions related to use of the NANOSOFT® include small risk of micro bleeding (self-limiting blood dotting, ceasing within normal clotting time). Other potential device related adverse events may rarely include local infection at the injection site.

STORAGE AND HANDLING

- NANOSOFT® should be stored at -10 to 30°C (14 to 86 °F).
- The NANOSOFT® is individually packed in a soft blister, sterilized by ethylene oxide (EtO).
- Do not use the NANOSOFT®, in the event that the soft blister is either open and/or damaged.
- Dispose of all devices in accordance with national/institutional guidelines, and in their absence, discard in a sharps container.




SUPPLY METHOD

The device is supplied as a single sterile NANOSOFT® device, packed in a soft blister for single use.

EN

NANOSOFT® MICRONEEDLES

	MANUFACTURER		KEEP DRY		MEDICAL DEVICE
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/ EUROPEAN UNION		DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE		SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE, STERILIZED USING EtO
	USE-BY/EXPIRATION DATE		DO NOT REUSE		PACKAGING UNIT
	LOT/BATCH NUMBER		CONSULT INSTRUCTIONS FOR USE		TEMPERATURE LIMIT -10°C / 14°F to 30°C / 86°F
	CATALOGUE NUMBER		NON PYROGENIC		PRESCRIPTION ONLY (USA)
	MODEL TYPE		SINGLE STERILE BARRIER SYSTEM, STERILIZED USING EtO		CE CONFORMITY MARKING 0483
	DATE OF MANUFACTURE		UNIQUE DEVICE IDENTIFIER		FRAGILE, HANDLE WITH CARE
	DO NOT RESTERILIZE		KEEP AWAY FROM SUNLIGHT		

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